

ESBL RESISTANCE DETECTION

NG-Test® CTX-M MULTI is a qualitative lateral flow immunoassay test that provides rapid detection of five main groups of CTX-M type enzymes of extended spectrum β -lactamases (ESBL) produced by Enterobacterales from a bacterial colony.

NG-Test®CTX-M simultaneously detects enzymes of CTX-M 1-2-8-9 and 25, including the most clinically relevant variants, from a bacterial colony in less than 15 minutes

NG-Test®CTX-M MULTI

